

Stability 4th



PART Name. :	BOBBIN WITH PIN	Instrument Name:	micrometer	Appr Name. :	ภณณ
PART NO. :	MK126EF447	Serial Number:	25648	Trials.	3
Characteristic/Specification	ยาว 0.0±2 mm	Resolution:	0.001		

LCL _X	-0.0194	CL _X	-0.0089	UCL _X	0.002	CL _R	0.1516	UCL _R	0.3903
------------------	---------	-----------------	---------	------------------	-------	-----------------	--------	------------------	--------

รอบ วันที่		1	2	3	4	5	6	7	8	9	10	11	12
		25-Jan	18-Feb	19-Mar	25-Apr	27-May	16-Jun	20-Jul	17-Aug	6-Sep	11-Oct	23-Nov	24-Dec
Reference value		0.348	0.584	0.350	-0.600	1.200	0.550	-0.854	-0.099	0.664			
Measured Values	1	0.290	0.556	0.290	-0.560	1.340	0.470	-0.800	0.020	0.590			
	2	0.410	0.554	0.410	-0.680	1.170	0.500	-0.920	-0.110	0.750			
	3	0.340	0.640	0.343	-0.580	1.270	0.640	-0.840	-0.210	0.660			
	Ave	0.347	0.583	0.348	-0.607	1.260	0.537	-0.853	-0.100	0.667			0.2423
	R	0.120	0.086	0.120	0.120	0.170	0.170	0.120	0.230	0.160			0.1440

รอบ		1	2	3	4	5	6	7	8	9	10	11	12
Bias	1	-0.058	-0.028	-0.06	0.04	0.14	-0.08	0.054	0.119	-0.074			
	2	0.062	-0.03	0.06	-0.08	-0.03	-0.05	-0.066	-0.011	0.086			
	3	-0.008	0.056	-0.007	0.02	0.07	0.09	0.014	-0.111	-0.004			
	Ave	-0.001	-0.001	-0.002	-0.007	0.060	-0.013	0.001	-0.001	0.003			0.0042
	R	0.120	0.086	0.120	0.120	0.170	0.170	0.120	0.230	0.160			0.1440

Result	Bias	Pass	Pass	Pass	Pass	Reject	Pass	Pass	Pass	Reject			
	Repeat	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass			



A2 and D3 and D4 are constant Value received from staitstical table				New	$\overline{R} =$	0.1440
Trials.	A2	D3	D4		$\overline{\overline{X}} =$	0.0042
2	1.880	0.000	3.267		$UCL_R = \overline{R} \times D_4 =$	0.3708
3	1.023	0.000	2.575		$LCL_R = \overline{R} \times D_3 =$	0.0000
					$UCL_x = \overline{\overline{X}} + A2 * R$	0.152
					$LCL_x = \overline{\overline{X}} - A2 * R$	-0.143

